

New
solution !

ITREMOR®

Drive safer in your wheelchair





ITREMOR®

Do you have tremors and difficulties with driving your wheelchair?

Do your tremors inhibit you to drive a wheelchair?

We have a solution for you!


iTremor® is a movement stabilizer which ensures that unwanted tremors on the joystick are not converted to unwanted movement of your wheelchair. It is inserted between the joystick and the wheelchair control.



The adjustment of your movements is easy via your Smartphone with iTremor® application. Subsequent adjustments are made via your medical facility or our field service. If your movements change over time, you can easily recalibrate it according to your preferences.

The iTremor® can be connected to different joysticks or chin control and different wheelchair models.

*Make your driving safer
and get a higher independence!*



iTremor® is already helping people with:

- **Parkinson's disease**
- **Muscular dystrophy**
- **Spasticity**
- **Muscular atrophy**
- **ALS**

iTremor® has also been designed for users with strong spasticity. It means that the solution has additional specific setting possibilities for speed reduction and shutdown options. In this case, iTremor® works best with the Heavy Duty Joystick and can be connected to different wheelchair models.

It helps patients with strong involuntary movements avoid dangerous situations and thus make driving safer.

Suitable joysticks



Produkt name	MicroPilot	Micro	Multi	Compact	Heavy Duty
Actuation force	1g	3.3g	8g	50/220g	750g
Actuation path	0.1mm	1mm	5mm	17mm	40mm

iTremor® is CE approved

If you are interested in iTremor®

- **Call us at +45 9635 4480**
- **Send us e-mail at info@tko-technology.dk**
- **Meet us at our office**



TKS Technology A/S
Niels Jernes Vej 10
DK9220 Aalborg
+45 9635 4480
info@tko-technology.dk
www.tko-technology.dk